





HIV-1/HIV-2 Antibody Test



- + Easy to use, early detection <sup>1</sup>
- + Clear and accurate results
- + Enhances patient flow
- + Expedites linkage to care





biolytical.com • insti.com













# INST!) HIV-1/HIV-2 Antibody Test

Early IgM detection <sup>1</sup>


**IgG Capture Procedural Control** 

No timers required

Lancet, pipette & alcohol swab included

"The ease of use of INSTI® and almost instant results make this test an accurate option for settings where recommended laboratory is not feasible... " - Adams, CDC Study 1

## Sample & pour



## **Product Information**

Information Type	Product Details
Method	Flow-through
Time to results	As little as 60 seconds
Storage conditions	35.6-86 °F (2-30 °C)
Shelf life	15 months
Sample type	Fingerstick, venous whole blood, plasma or serum


## Ordering Information

#### Product

INSTI® HIV-1/HIV-2 Antibody Test

INSTI® HIV-1 Positive/HIV-2 Positive HIV Negative Controls (3 x 1.0 mL vials)

## Read results immediately!



## Accuracy <sup>2</sup>

Sample Type	Sensitivity	Specificity
Venous Whole Blood	99.4%	99.7%
Plasma	99.5%	99.96%
Fingerstick Blood	99.6%	99.3%
Serum	99%	100%



## **Trusted Technology That Works.**

Confidence. The protein-A control feature serves as a built-in internal control and confirms correct sample addition, test procedure and proper test performance.

### Distributed by:

#### **Una Health Ltd**

T: 01782 575180

E: enquiries@unahealth.co.uk

W: www.unahealth.co.uk



#### References

<sup>1</sup> Adams S, et al. Performance evaluation of the point-of-care INSTI® HIV-1/2 Antibody Test in early and established HIV infections. Journal of Clinical Virology (2017) 91: 90-94.

<sup>2</sup> INSTI HIV-1/HIV-2 Antibody Test Package Insert (CAN), 51-1028 (J).