

DIAGNOSE WITH CONFIDENCE

IMMUVIEW[®]

***S. pneumoniae
and Legionella***

Identify both with
just one test



For urine and CSF samples



Una Health Ltd
Unit 6 & 7 Fitzgerald Way
Scotia Road Business Park
Stoke-on-Trent, ST6 4HN
T +44 (0)1782 575180

unahealth.co.uk

una
health

IMMUVIEW® LIMIT OF DETECTION


Limit of detection ImmuView® and BinaxNOW™, 1

Dilution

Pneumococcal antigen (mg/ml)	ImmuView®	BinaxNOW™
1×10^{-9}	Positive	Positive
1×10^{-10}	Positive	Negative
1×10^{-11}	Positive	Negative
1×10^{-12}	Negative	Negative

"ImmuView® is more sensitive than BinaxNOW™ for the detection of purified pneumococcal antigen"¹

IT ONLY TAKES 3 SIMPLE STEPS TO GET RESULTS



Why ImmuView®?


2 results in 1 test

Similar symptoms
– different treatment

Get results in just 15 minutes

1. Adapted from, Lindsay et. All, Evaluation of the rapid lateral flow ImmuView® *S. pneumoniae* and *L. pneumophila* urinary antigen test (UAT) for the detection, poster presented at ESGLI 2017

Product details
ImmuView® *S. pneumoniae* and *L. pneumophila* Urinary Antigen Test, Article No. 95389 . Kit Size: 22 tests.



DIAGNOSE WITH
CONFIDENCE