

# Point of Care CBC analyzer

- Single finger prick.
- Two drops of blood.
- Accurate results in minutes, not days.



🕒 45 Seconds

## 1. Sample Preparation

Using OLO's single-use test kit, add 2 drops of capillary or venous blood to the disposable cartridge sample



🕒 10 Minutes

## 2. Automated Scanning

Insert cartridge and wait for machine-vision scanning to be performed (blood sample never comes into contact with the device)



🕒 Immediate

## 3. Results Visualization

Visualize all results on OLO's touchscreen and transfer data to the patient record



Factory calibrated for easy set-up

Remote installation, training and servicing

Capillary and venous blood sample

No maintenance required

No external reagents needed

No liquid waste


Uniquely designed for unprecedented times.  
And, it just works.

Intended for audience outside of the United States | CMR0133 rev 1.0  
All rights reserved | ©2020 Copyright Sight Diagnostics®

sight  
[www.sightdx.com](http://www.sightdx.com)

Compact, Convenient.  
Cost-effective.

**Sight OLO<sup>®</sup> has been placed  
in an isolated environment  
designated for supporting  
COVID-19 patients**



**Size Specification**

Height: 323mm

Width: 284mm

Depth: 254mm



OLO's compact design makes it easy to set-up quickly in a setting reserved for testing only contagious blood samples. Eliminating the risks of processing these samples alongside regular lab samples is uniquely advantageous in outbreaks of highly infectious diseases.

**Dr Ram Doolman**

Director of Laboratories Division and Automated  
Mega-Laboratory at Sheba Medical Centre at Tel HaShomer



**Contact Sight Diagnostics**

**Contact Sight Diagnostics**  
[www.sightdx.com](http://www.sightdx.com)

**Email**  
[sales-ww@sightdx.com](mailto:sales-ww@sightdx.com)

CE Marked for Point of Care. TGA Approved. FDA 510 (k) cleared for moderately complex settings.  
For complete indications for use and safety information of Sight OLO<sup>®</sup>, please visit [www.sightdx.com](http://www.sightdx.com).

**Once you've seen how it works,  
you'll wonder how you ever worked without it.**

Intended for audience outside of the United States | CMR0133 rev 1.0  
All rights reserved | ©2020 Copyright Sight Diagnostics<sup>®</sup>

sight  
[www.sightdx.com](http://www.sightdx.com)