





Point of care HbA1c testing has shown to improve diabetes management, patient's commitment to treatment and patient satisfaction. These are leading to improved glycemic control and reduced risk of harmful and expensive comorbidities. It also has the potential to reduce the healthcare costs.1-3



### Easy sample collection and test procedure

- One-step sample collector
- Results in less than 6 minutes
- 5 minutes sample stability in the cuvette
- Automated & portable QuikRead go instrument



#### Reliable HbA1c results at the point of care

- Certified Silver by the IFCC
- Multiple **EQA schemes** available
- HbA1c results available during patient consultation



# Flexibility in Testing

- Suitable for **all patients**, even in paediatrics
- Small finger prick sampling (1 µl)



# **Full Portability & Connectivity**

- Light-weight device with carrier bag
- Can be used on **battery** or mains power
- **High battery performance** 50/100 tests
- Extended **connectivity** options (HIS/LIS)

#### Did you know?


The QuikRead go HbA1c method is certified Silver by

IFCC certification is a requirement of NICE guidelines.

# Efficacy and simplicity with point of care HbA1c testing

Enhancing diabetes management pathways | Streamlining workflows | Improving patient's commitment to treatment and monitoring

**HbAIc testing at the point of care** has the potential to **reduce costs** by **streamlining the workflow** in healthcare and by **lowering the risk** of harmful and expensive **complications of diabetes**. Point of care HbAIc has shown to **improve** patient satisfaction and **patient's acceptance and commitment to treatment** and **monitoring**. The risk of harmful complications decreases with the **improved glycemic control**. The point of care result is ready during the appointment and no separate visits are needed. Treatment plan can be discussed immediately, and no time is lost. <sup>1-5</sup>





#### References

- 1. Miller CD, Barnes CS, Philips LS, et al. Rapid Alc availability improves clinical decision-making in an urban primary care clinic. Diabetes Care 2003; 26(4):1158-1163.
- 2. Schnell O, Crocker JB, Weng J. Impact of HbA1c Testing at Point of Care on Diabetes Management. J Diabetes Sci Technol 2017; 11(3):611-617
- 3. American Diabetes Association. Diabetes Care 2020; 43(Supplement 1): S66-S76.
- 4. Patzer KH, Ardjomand P, Gohring K, et al. Implementation of HbAlc Point of Care Testing in 3 German Medical Practices: Impact on Workflow and Physician, Staff and Patient satisfaction. J Diabetes Sci Technol. 2018;12(3):687-694
- 5. Rosa LS, Mistro S, Oliveira MG et al. Cost-effectiveness of Point-of-Care Tests in Primary Care Setting. Front. Pharmacol. 2021 Jan 19;11:588309
- 6. NICE, Type 2 diabetes in adults: management, NICE Guideline [NG28], https://www.nice.org.uk/guidance/ng28/chapter/recommendations