

ImmuView[®] Reader 2.0

Automated reading and
result interpretation

The ImmuView[®] Reader 2.0 is a compact, easy-to-use instrument for automated reading and result interpretation.



It is compatible with the *C. DIFF QUIK CHEK COMPLETE[®]* assay to provide:

- **Objectivity:** Precise interpretation of results
- **Connectivity:** LIS integration capability
- **Traceability:** On-board data storage
- **Simplicity:** User-friendly touchscreen interface
- **Small footprint:** 124 x 118 x 142 mm (WxHxD)



ImmuView® Reader 2.0

Automated reading and result interpretation

Benefits

- **Ethernet port** for network connection to LIS.
- **Simple user interface** - LCD touchscreen is compatible with gloves and provides visibility in a variety of lighting conditions. The graphic user interface is easy to navigate with an intuitive menu structure.
- **Barcode Scanner and Printer** - easy plug and play connectivity for USB printer and barcode scanner.
- **Straightforward workflow** - place the test in the cartridge drawer, close the drawer, and press the test button.
- **Store up to 999 tests** - you can always retrieve previous results. The ImmuView® Reader will inform you before it reaches the storage limit.
- **Integrated LIS/LIMS compatibility** with automatic data transfer for prompt and accurate reporting of test results.

