

## Lepzi® PLGF Test

Aiding diagnosis and prognosis of pre-eclampsia



**Point of care PLGF testing  
direct from whole blood**

The Lepzi PLGF test is a next generation point of care test that can be performed without the need for laboratory facilities, in just 15 minutes. No centrifuging required.



Aids diagnosis and  
prognosis of  
pre-eclampsia  
between 20- and  
37-weeks' gestation



**Fast results in  
15 minutes**




**UKCA marked**



**Comparable  
performance to NICE  
recommended tests**

### Test procedure



**Step 1:** 200µl whole  
blood added to cartridge



**Step 2:** Cartridge  
inserted into Lepzi  
Reader



**Step 3:** Read results  
on-screen, or press  
'print' for a hard copy

### Product information

Information	Product details
Method	Fluorescence immunoassay
Storage	2 - 8°C
Shelf life	12 months
Sample type	EDTA-anticoagulated blood or plasma